

# Refreshing energy!!!

It helps change your scalp and skin far clearer and brighter and increase your skin resilience by improving skin moisturizing and shine.

## Beading Action & Adipo Therapy System

### Mu:Beader™ -Bp

It carries out from cleansing, exfoliation of aged keratin up to removal of sebum painlessly using micro-bubbles generated by a new concept of natural spa skin care and a high-frequency massager and has multiple functions to improve skin moisturizing and luster and massages the skin simultaneously using a high-frequency massager.



### Special Tip

### The characteristics

Removal of horny substances and supply of moisture at the same time.

Mu:Beader™-Bp is designed to promote exfoliation with fine bubbles to remove aged dead skin cells and sebum and, at the same time, inject highly concentrated, nutritious serum so as to make skin and scalp moist, smooth, and lustrous.



The special tip of the Mu:Beader™Bp consists of two suction holes and one serum discharge hole. Its specially-designed spiral tip is intended to hold a product long enough to soften dead cells of the scalp and skin, provide nutrients, and absorb them in a vacuum state to maximize its peeling and extrusion effects.

Safety : neither downtime nor side-effect

Because treatment does not cause any damage to skin and scalp, patients can lead a normal life right after treatment. There are no side effects such as injuries, pain, burns, rashes, etc.

Multiple effects on skin care

Makeup absorption rate and skin moisturizing and luster will be enhanced by the action of biological heat generated by high-frequency.

### A hand-piece of high-frequency



There are two types of hand-piece; one is 25mm and the other 40mm, which are all using high-frequency and used selectively depending on the size of treatment areas.

## Feature of high- frequency

- Decrease of Subcutaneous Fat: Body Shape Correction
- Lively electronic activity is made around handpiece that comes into contact with skin. So it's effective for increasing the temperature of subcutaneous fat which is close to surface.

It promotes the circulation of the blood and the regeneration of skin, and melts subcutaneous fat cells. It's also effective for body-lifting.

## User Manual for Mu:Beader™ -Bp

### 01 Cleansing

Wipe wastes 2~3 times off the surface of the face with a cleansing cream.

### 02 Treatment(removal of keratin and sebum)

First spread the selected fine bubble over the face and then divide the face into 3~5 areas to start treatment for about 2 minutes.

### 03 Treatment of high-frequency

Simple to use as it has no plate. Make treatment as if you give a massage on the desired area after applying a designated cream excellent in current carrying capacity to the area to be cared.



### 04 Washing

Wash and wipe out the face cleanly with a spongy after treatment.

## Cautions in the treatment by Mu:Beader™ -Bp

### 01. Hand-pieces of vacuum compression type

Treatment shall be provided in the only condition that a handpiece is in fully contact with the skin of a patient. Otherwise, the discharge and suction of bubbles will not be made normally.

### 02. Temporary redness phenomenon

The redness in the skin occurred by the skin irritation during the treatment is normal and temporary and returns to original skin tone in approximately 10 ~ 30 minutes.

### 03. Treatment order of neck – chin – cheeks – forehead – nose

Provide treatment as if the entire face is lifted from the lower part to the upper part along the skin texture.

### 04. A hand-piece of high-frequency massage

Simple to use as it has no plate. Make treatment as if you give a massage on the desired area after applying a designated cream excellent in current carrying capacity to the area to be cared.